

LC/MS and ELSD Analytical Procedures for Ionic Liquids

Ted Chang, Michael Piquette, Liying Du, Cytec Industries Inc., Stamford, CT/USA

D. Nucciarone, I. Water-Protas, Cytec Canada Inc., Niagara Falls, Ontario/CDN

Introduction

Phosphonium ionic liquids (IL) are an important class of ionic liquids with enhanced stability with respect to degradation under various conditions and thermal stability.

In this report, we present an LC/MS procedure for the analysis and purity determination of mixtures of phosphonium IL. The method reported is a robust procedure providing excellent peak resolution and high detection sensitivity.

Since most cationic species of phosphonium type ionic IL are not sensitive to UV detection, this procedure utilizes two mass spec ionization techniques (Electrospray and APCI). A third HPLC detection technique (ELSD) is also evaluated.

Design of Experiment


1. Selection of representative IL: four phosphonium IL (two low MW and two high MW) and two imidazolium.
2. Evaluation of the separation efficiency of HPLC column (column type, mobile phase, gradient conditions).
3. Selection of optimum HPLC conditions.
4. Comparison of LC/MS and LC/ELSD procedures: resolution and sensitivity.
5. Identification of anionic species from dimerized MS ion.
6. Analysis of residual chloride to 10ppm level.

Conclusions


1. An LC/MS procedure is established for phosphonium type IL.
2. This procedure is also applicable for imidazolium type IL.
3. Sensitivity is 1 ppm (0.1% using 0.1% sample solutions).
4. LC/ELSD gives better peak resolution.
5. No LC/MS procedure specific for anionic species has been developed, however, anions can be identified from dimerized cluster ions.
6. IC can determine residual chloride in non-chloride IL to 10 ppm using 1% sample solutions.

LC/MS and ELSD Procedures

LC/MS-ESI and HPLC-ELSD Chromatograms



Mass Spectra of IL Cationic Species



Observed Cations

- P1 Triisobutyl(methyl)phosphonium, m/z = 217
- P2 Tetrabutylphosphonium, m/z = 259
- P3 Tributyl(tetradecyl)phosphonium, m/z = 399
- P4 Trihexyl(tetradecyl)phosphonium, m/z = 483
- Im1 1-Ethyl-3-methyl-1H-imidazolium, m/z = 111
- Im2 3-Butyl-1-methyl-1H-imidazolium, m/z = 139

Analytical Procedure

Mass Spectrometer: Applied Biosystems Mariner LC/MS system
Mass Spec ionization techniques: Electrospray and APCI

HPLC: Agilent HP 1100
HPLC detector: ELSD


HPLC Conditions:
Column: Varian MetaSil Basic, 250 x 4.6mm
Aqueous mobile phase: 1% HAc, 0.2% NH₄Ac, pH = 3.9
1% HAc only, pH = 2.8
Organic mobile phase: Methanol
Flow rate: 0.5 mL/min.
Gradient: 50-100% MeOH (20min.) + 15 min.
Column temperature: 30°C
Sample size: 5 µL of 0.01% in MeOH

Relative Signal Response from Different Detectors (normalized to P2)


Detector	P1	P2	P3	P4
APCI	0.42	1	1.57	1.67
Electrospray	0.85	1	0.86	0.67
ELSD	0.32	1	3.00	1.70

Evaluation Criteria


Selection of Column




Effect of pH of Aqueous Phase



Effect of Gradient System




Linearity of LC/MS



P1: Green
P2: Blue

Linearity of ELSD



P3: Black
P4: Red


Results of Linearity Study

1. LC/MS has excellent sensitivity and linearity.
 - Sensitivity is 1 ppm.
 - Can detect 0.1% impurity with a 0.1% sample solution.
 - Linearity is 2 orders of magnitude.
2. ELSD gives better peak resolution.

Identification of Anions

MS CLUSTER IONS

Identification of Anionic Species using Dimerized Cluster Ion Three IL with the Same Cation, Different Anions




Identification of Anionic Species

- All three IL have the same cation (m/z = 483)
- Anion of P6 is confirmed as m/z = 289
- Anion of P7 is confirmed as m/z = 280
- Anion of P8 is confirmed as m/z = 66
- Impurity in P6 is chloride, m/z = 35

Ion Chromatography Linearity for Chloride

(Sensitivity: 0.1ppm = 10ppm in sample using 1% solution)



Residual Chloride in non-Chloride Ionic Liquids (using 1% solution)

IL abbr.	Residual Cl in sample
P1	80 ppm
P5	50
P6	160
P8	630
P9	890
P10	280
P11	540
P12	170
Im3	10
Py1	500


List of Analyzed Ionic Liquids

Ionic Liquids used in HPLC Study (LC/MS and ELSD)

- P1 CYPHOS IL® 106 Triisobutyl(methyl)phosphonium tosylate
- P2 CYPHOS IL 164 Tetrabutylphosphonium chloride
- P3 CYPHOS IL 167 Tributyl(tetradecyl)phosphonium chloride
- P4 CYPHOS IL 101 Trihexyl(tetradecyl)phosphonium chloride
- Im1 1-Ethyl-3-methyl-1H-imidazolium chloride
- Im2 3-Butyl-1-methyl-1H-imidazolium chloride

Ionic Liquids used for Anion Study (MS and IC)

- P5 CYPHOS IL® 163 Tetrabutylphosphonium bromide
- P6 CYPHOS IL 104 Trihexyl(tetradecyl)phosphonium (ditrimethylpentyl)phosphinate
- P7 CYPHOS IL 109 Trihexyl(tetradecyl)phosphonium bis(trifluoromethane)sulfonylamide
- P8 CYPHOS IL 105 Trihexyl(tetradecyl)phosphonium dicyanamide
- P9 CYPHOS IL 162 Tributyl(hexadecyl)phosphonium bromide
- P10 CYPHOS IL 107 Tetrabutylphosphonium dibutylphosphate
- P11 CYPHOS IL 102 Trihexyl(tetradecyl)phosphonium bromide
- P12 CYPHOS IL 110 Trihexyl(tetradecyl)phosphonium hexafluorophosphate
- Im3 1-Ethyl-3-methyl-1H-imidazoliumtrifluoromethanesulfonate
- Py1 Pyridinium p-toluenesulfonate



References

1. C. Bradaric, A. Downard, C. Kennedy, A. Robertson and Y. Zhou, *Green Chemistry* 5, 143 2003
2. T.T. Chang and M.J. Piquette, *Dechema Conference on Green Solvents*, Bruchsal, Germany, 2002.
3. P. Stepnowski, A. Muller, P. Behrend, J. Ranke, J. Hoffmann and B. Jastorff, *J. Chromatog. A*, 993, 2003.
4. R. Del Sesto, A. Robertson, et. al., *J. Organometallic Chem.*, 690 (10), 2005

Acknowledgements

Eduardo Kamenetzky Cytec Industries Inc., Stamford, CT USA
Al Robertson Cytec Canada Inc., Niagara Falls, Ontario, Canada